

# GC Initial™ IQ SQIN Powders



This chart is intended as a reference guide only. The "Form & Texture" Liquid is a special and dedicated liquid to be used only with the SQIN powders.

# GC Initial™ IQ SQIN Powders



The SQIN ceramics are a refined mixture of feldspar-based glasses dedicated to the Initial IQ ONE SQIN Concept – a paintable color and form ceramic technique for the aesthetic realization of monolithic or minimal buccally reduced dental frameworks made of zirconium dioxide or lithium disilicate e.g. GC Initial LiSi Press and GC Initial Zirconia Disk. For more tips and tricks please consult the dedicated GC Initial IQ ONE SQIN Technical Manual.

V-Shade		A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
Dentin Body Shades	4	Body A					Body B				Body C				Body D			
Enamel	4	E58	E58	E59	E59	E60	E57	E59	E59	E59	E60	E59	E59	E60	E60	E59	E59	
Transpa Opal Booster	1	TO Booster																
Bleach	2	BL-D (Bleach Dentin)									BL-E (Bleach Enamel)							

